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Instruction Manual

The Domotherm® OT Digital Infrared Ear Thermometer is a high quality product incorporating the latest technology and tested in accordance with international standards. With its unique technology, the Domotherm® OT can provide a stable, heat-interference-free reading with each measurement. The instrument performs a self-test every time it is switched on to always guarantee the specified accuracy of measurements.

Why measure in the ear?

The goal of thermometry is to measure core body temperature which is the temperature of the vital organs. Ear temperatures accurately reflect core body temperature, since the eardrum shares blood supply with the temperature control center in the brain, the hypothalamus. Therefore, changes in body temperature are reflected sooner and more accurately in the ear than at other sites.

As with any thermometer proper technique is critical to obtaining accurate temperatures. Therefore, read the use instructions carefully and thoroughly.

UEBE Medical GmbH

Table of Contents

1. The Advantages of your Domotherm® OT Digital Infrared Ear Thermometer
2. Important Safety Instructions
3. Product Description
4. How the Domotherm® OT Digital Infrared Ear Thermometer Measures Ear Temperature
5. Control Displays and Symbols
6. How to Reload a New Probe Cover
7. Direction for Use
8. Changing from Fahrenheit to Celsius and vice-versa
9. Cleaning and Storage
10. Error Messages
11. Replacing the Battery
12. Technical Specifications
13. Guarantee
14. www.uebe.com

1. The Advantages of your Domotherm® OT Digital Infrared Ear Thermometer

Wide Range Measurement

Additionally to the requested range of 35,5 °C bis 42,0 °C (95,9 °F bis 107,6 °F) for body temperature measurement Domotherm® OT offers you a wide range measurement of 0 °C to 100 °C (32 °F to 212 °F).

Measurement in 1 second

The innovative infrared technology allows measurement of ear temperature in only 1 second.

Accurate and reliable

Due to the unique probe assembly construction, the advanced infrared sensor, and the complete calibration process this unit can offer a very accurate and reliable ear temperature measurement.

Gentle and Easy to Use

- Special ergonomic design enables simple and easy use of the thermometer.
- The Domotherm® OT Digital Infrared Ear Thermometer can be used without interference to daily lifestyle. A measurement can be taken even while a child is sleeping.
- The Domotherm® OT Digital Infrared Ear Thermometer is more pleasant to use for children than a rectal thermometer or an oral thermometer.

Auto-Display Memory

The product displays the last reading automatically for 2 seconds when the unit is switched ON.

Safe and Hygienic

- No risk of broken glass or mercury ingestion.
- Completely safe for use on children.
- Disposable probe covers make Domotherm® OT completely hygienic.

Fever Alarm

10 short beeps alert the patient that he/she may have fever.

NOTE:

When consulting your physician, communicate that the Domotherm® temperature is a temperature measured in the ear and if possible, note the individual's normal Domotherm® temperature range as additional reference.

2. Important Safety Instructions


- Never use the thermometer for purposes other than those it has been intended for. The ear thermometer is only intended for temperature measurements in the auditory canal. It should not be used for rectal, oral or underarm temperature measurements. Please follow the general safety precautions when using on children.
- Always use the thermometer with a new undamaged Probe Cover for each measurement to prevent infection. Only Domotherm® OT branded probe covers can ensure that you get an accurate measurement from the Domotherm® OT.
- Never immerse the Domotherm® OT Digital Infrared Thermometer into water or other liquids (not waterproof). For cleaning and disinfecting please follow the instructions in the „Cleaning and Storage“ section.
- Keep the instrument and the probe covers away from direct exposure to the sun and keep it in a dustfree, dry area at the temperature between 10 °C - 40 °C (50 °F - 104 °F).
- Do not use the thermometer if there are signs of damage on the measuring tip or on the instrument itself. If damaged, do not attempt to repair the instrument! covers from Uebe (Art.no. 8601, PZN 4084666)
- Earwax in ear canal may cause a lower temperature reading. Make sure subject's ear canal is clean to ensure an accurate reading.
- This Domotherm® OT Digital Infrared Ear Thermometer consists of high-quality precision parts. Do not drop the instrument! Protect it from severe impact and shock. Do not twist the instrument and the measuring probe!

WARNING:

- **Please keep the probe covers out of the reach of children.**
- **Use of this IR thermometer is not intended as a substitute for consultation with your physician.**
- **Thermometer is not waterproof! Please NEVER immerse into liquids!**
- **Domotherm® OT is not appropriate for measuring waking temperature.**

3. Product Description

- ① Probe Cover
- ② Probe
- ③ Start button
- ④ LCD Display
- ⑤ O/I button
- ⑥ Battery Cover

4. How the Domotherm® OT Digital Infrared Ear Thermometer Measures Ear Temperature

The Domotherm® Digital Infrared Ear Thermometer measures infrared energy radiated from the eardrum and the surrounding tissue. This energy is collected through the lens and converted to a temperature value. The measured reading obtained directly from the eardrum (Tympanic Membrane) can ensure the most accurate ear temperature. Measurements taken from the surrounding tissue of the ear canal generate lower readings and may result in misdiagnosis of a fever.

To avoid an inaccurate measurement:

- First slip on a probe cover.
- Switch on the thermometer by pressing the O/I button.
- After one beep is heard (and the temperature scale icon is flashing), straighten the ear canal by gently pulling the middle of the ear back and up.
- Place the probe firmly into the ear canal, press the Start button and keep the probe in the ear until the thermometer beeps to identify the completion of the measurement.

The Domotherm® Digital Infrared Ear Thermometer has been clinically tested and proven to be safe and accurate when used in accordance with its operating instruction manual.

5. Control Displays and Symbols

LCD-Display	Display Meaning	Description
	All segments displayed	Press the O/I button to run on the unit, all segments will be shown for 2 seconds.
	Memory	The last reading will be shown on the display automatically for 2 seconds.
	Ready	The unit is ready for the measurement, the °C or °F icon will keep flashing.
	Measurement complete	The reading will be shown on the LCD display with the °C or °F icon flashing, the unit is ready again for the next measurement.
	Low battery indication	When the unit is turned on, the battery icon will keep flashing to remind the user to replace the batteries.

6. How to Reload a New Probe Cover

(1) Place a probe cover onto the storage case hole with paper side upwards.

(2) Take the unit, vertically penetrate the probe into the center part of the probe cover.

(3) Completely push the probe into the probe cover holder hole.

(4) After feeling a slight click, take out the probe with cover attached tightly.

NOTE:

- In order to avoid cross-contamination, please reload a new probe cover for each measurement.
- Please check if the probe cover is fitted on firmly before use; If the probe cover is broken, discard the probe cover and reload a new one immediately.

7. Directions for Use**Important:**

Prior to every measurement, fit a new undamaged Probe Cover on the measuring probe. Without a probe cover the measuring probe will be damaged (no warranty claim). Failure to do so may result in incorrect temperature measurement.

1. Press the O/I button. The LCD is activated to show all segments for 2 seconds.
2. The last measurement reading will be shown on the display automatically for two seconds with the „M“ icon.
3. When the °C or °F icon is shown, a beep sound is heard and the thermometer is ready for the measurement.
4. Straighten the ear canal by pulling the ear up and back to give a clear view of the eardrum.
 - For children under 1 year: Pull the ear straight back.
 - Children aged 1 year to adult: Pull the ear up and diagonally back.
5. While tugging the ear, insert the probe snugly into the ear canal and press the „START“ button. Release it, when you hear a beep sound. This is the reminding signal that confirms the end of measurement.
6. Remove the thermometer from the ear canal. The LCD displays the measured temperature.

NOTE:

10 short beeps will sound when the temperature is higher than 37.5 °C in order to alert the patient that he/she may have fever.

7. Replace the probe cover after each measurement. To do this, please follow the instruction in point 6 „How to reload a new and clean probe cover?“.
8. In order to assure the accurate readings, please wait at least 30 seconds after every measurement.

NOTE:

- For an infant, it is best to have the child laying flat with his head side-ways so the ear is facing upwards. For an older child or adult, it is best to stand behind and slightly to the side of the patient.
- Always take the temperature in the same ear, since the temperature readings may be different from the right ear and left ear.
- Please wait for a few minutes to take the ear temperature after sleeping.
- In the following situations it is recommended that three temperatures in the same ear be taken and the highest one taken as the reading:
 - 1) New born infants in the first 100 days.
 - 2) Children under three years of age with a compromised immune system and for whom the presence or absence of fever is critical.
 - 3) When the user is learning how to use the IR thermometer for the first time until he/she has familiarized himself/herself with the instrument and obtains consistent readings.

8. Changing from Fahrenheit to Celsius and vice-versa

The Domotherm® OT Digital Infrared Thermometer can display temperature measurements in either Fahrenheit or Celsius. To switch the display between °C and °F, simply turn OFF the unit, press and hold the Start button for 5 seconds; When you stop pressing the Start button after 5 seconds, the current measurement unit (°C or °F icon) will be flashed on the display. Toggle the measurement unit between °C and °F by pressing the start button again. When the measurement unit has been chosen, wait for 5 seconds and the unit will enter the ready for measuring mode automatically.

9. Cleaning and Storage

Use an alcohol swab or cotton swab moistened with alcohol (70 % Isopropyl) to clean the thermometer casing and the measuring probe. Ensure that no liquid enters the interior of the thermometer. Never use abrasive cleaning agents, thinners or benzene for cleaning and **never immerse the instrument in water or other cleaning liquids**. Take care not to scratch the surface of the LCD.

Remove the battery from the instrument if it is not required for extended periods of time in order to avoid damage to the thermometer resulting from a leaking battery.

10. Error Messages

Display/Problem	Display Meaning	Possible cause and fault remedy
	Measured Temperature too high	Displays „H“ when measured temperature higher than 100.0 °C or 212 °F.
	Measured temperature too low	Displays „L“ when measured temperature lower than 0 °C or 32.0 °F.
	Ambient temperature too high	Displays „H“ in conjunction with the „▲“ when ambient temperature is higher than 40.0 °C or 104.0 °F.
	Ambient temperature too low	Display „L“ in conjunction with the „▼“ when ambient temperature is lower than 10.0 °C or 60.8 °F.
	Error function display	When system has malfunction.
	Blank display	Please check if the battery has been loaded correctly. Also check polarity (<+> and <->) of batteries.
	Dead battery indication	If the steady battery icon is the only symbol shown on the display, the batteries should be replaced immediately.

11. Replacing the Battery

The Domotherm® OT Digital Infrared Thermometer is supplied with one lithium battery, type CR2032. Replace with a new CR2032 battery when the flashing battery symbol appears on the LCD display. Using a screwdriver to loosen the screw from battery cover as shown, remove the battery cover and replace CR2032 battery.

Caution:

Batteries and technical devices do not belong in the domestic waste and must be deposited in the respective collection disposal point.

12. Technical Specifications

Type	Digital Infrared Thermometer
Measuring Range	0 °C to 100.0 °C (32.0 °F to 212.0 °F)
Accuracy	Laboratory: : ± 0.2 °C, 32.0 – 42.2 °C (± 0.4 °F, 89.6 – 108.0 °F) ± 0.3 °C, 0 ~ 31.9°C, 42.3 ~ 100.0 °C (± 2 °F, 32.0 ~ 89.5 °F, 108.1 ~ 212.0 °F)
Display	Liquid Cristal Display with indicating unit 0.1 °C (0.1 °F)
Acoustic	a. The unit is turned ON and ready for the measurement: 1 short „bi“ sound. b. Complete the measurement: 1 long beep sound. c. System error or malfunction: 3 short „bi“ sounds.
Memory	Auto-Display the last measured temperature
Nite Glow	a. The display will be lighted for 4 seconds when the unit is turned ON. b. The display will be lighted again for 5 seconds when the measurement has been completed.
Operating temperature	10 °C to 40 °C (50.0 °F to 104 F)
Storage/transport temperature	-25 °C to +55 °C (-13 °F to 131 °F)
Automatic Switch-off	Approx. 1 minute after last measurement has been taken.
Battery	CR2032 BATTERY (X1) at least 1000 measurements
Dimensions	141mm (L) x 26 mm (W) x 20 mm (H)
Weight	45g (with battery), 40g (w/o battery)

This device conforms to the requirements of the following standards:

EN 12470-5:2003

Clinical Thermometers – Part 5: Performance of infra-red ear thermometers (with maximum device)

EN 60601-1:1990

Medical electrical equipment – Part 1: General requirements for
+ A1:1991 safety
+ A2:1995

EN 60601-1-2:2001

Medical electrical equipment – Part 1-2: General requirements for safety –
Collateral standard: Electromagnetic compatibility – requirements and tests

According to the Medical Product User Act a biennial technical inspection is prescribed for professional users.
Please observe the applicable disposal regulations.

This device conforms to Directive 93/42/EWG of June 14, 1993 on Medical Devices und bears the CE mark "CE 0123" (TÜV SÜD Product Service).

Protection grade against electric shock: Type BF

13. Manufacturer's Guarantee

The equipment was manufactured and tested with great care. However, in the unlikely event of a defect being detected after delivery, we are prepared to give a guarantee in accordance with the following conditions:

1. During the guarantee period of 2 years from the date of purchase, we reserve the right to either repair any such defect at our expense or supply a perfect replacement unit (in any case the defective equipment must be returned to our factory).
2. Excluded from the guarantee are parts subject to normal wear and tear as well as damage caused by non-compliance with the operating instructions, improper handling (for example: unsuitable supply mains, breakages, leaking batteries) and/or dismantling of the unit by the purchaser. Furthermore, no claims for damages against us can be justified through the guarantee.
3. **Claims under the guarantee can only be made during the guarantee period by presentation of the purchase receipt.** In the case of claim under the guarantee, the equipment together with the purchase receipt and a written justification of the claim are to be sent to

UEBE Medical GmbH
Zum Ottersberg 9
97877 Wertheim
GERMANY

4. In the case of defectiveness of the goods, the contractual rights of the purchaser to claim against the supplier/salesman in accordance with § 437 BGB (The German Civil Code) are not limited by the guarantee.

14. www.uebe.com

Detailed user information about our Domotherm® thermometers, visomat® and visocor blood pressure monitors, CycloTest products as well as services can be found at www.uebe.com.

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© Copyright UEBE Medical GmbH