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GB

Congratulations for buying this high quality product. To use all advantages of this thermometer you should familiarize with the manual prior to first measurement. For later use keep the manual stored in a safe place. If you follow the instruction manual you will be happy with your digital fever thermometer for years.

ADVANTAGES

- **Maximum-Thermometer.** The highest measured value remains visible until device switches off.
- **Acoustic signal** indicates end of measurement (10 long beeps)
- **Measured value** after approx. 60-90 sec., depending on method.
- **Fever alarm** from 37.80 °C on (10 x 3 short beeps)
- **Extra large digital display** with two digits after decimal point.
- **Waterproof.** Complete cleaning bath possible.
- **Storage.** Last measurement memory for comparison with next measurement.
- **Automatic switch off** for battery protection (after approx. 10 minutes).
- **Longlife-Battery** Replacement battery available through good specialized dealers.

IMPORTANT INFORMATION

- Do not drop the device. It is not shock-proof.
- Before measurement check intactness of tip and display.
- This thermometer is determined only for measuring human body temperature.
- This thermometer contains small parts which can be swallowed by children. Don't leave children with this device uncontrolled.
- Do not leave the device in direct sunlight or near heat sources.
- Do not bend the thermometer.
- Do not open the device (except the battery cover).
- For cleaning use only water or special substances recommended in chapter „Cleaning and Disinfection“.

INSTRUCTIONS FOR USE

Press the on/off-button. After an acoustic signal, all symbols show on the display for a few seconds as a means of checking that the thermometer is functioning correctly. If available the last correctly measured temperature is displayed indicated by the „M“ symbol in the bottom right corner. Subsequently the device indicates an internal test value of 37.00 °C. After this the device is ready for measurement.

During the measurement the current temperature is shown indicated by a flashing „C“. When the acoustic signal sounds and the „C“ stops flashing then the temperature increase is less than 0.1 °C in between 16 seconds. This means the nearly final temperature if you use the oral or rectal measurement method. The temperature then changes only minimal. With the axillary method you have to keep the specified time of measurement and the acoustic signal can then sound several times.

The measured value remains visible for approx. 10 minutes. After that the device switches off itself automatically. To save battery energy switch off the thermometer after use.

The previously measured temperature value remains stored in any case until the next measurement.

MEASUREMENT METHODS

Oral measurement (in the mouth)
Be sure to position the sensor tip under the tongue. This measurement takes approx. 1 minute.

Rectal measurement (in the anus)
Gently insert the sensor tip 1 cm into the anus. This measurement takes approx. 1 minute.

Axillary measurement (under the arm)
Insert the sensor tip under the arm pit so that the sensor tip is in close contact with the skin. This measurement takes approx. 5-10 minutes.
From the medial point of view this method always provides inaccurate values. And so it is unsuitable for accurate measurement.

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GB

WHAT SHOULD YOU DO IF YOU HAVE A TEMPERATURE?

A normal body temperature is 36.9°C ±0.5°C. Temperatures from 37 to 38.5°C are called a moderate temperature and above 39°C a high temperature (rectal method). Natural methods are very suitable for lowering the temperature, for example:

- 1. Sponging with cold water (water 10°C lower than temperature measured).
- 2. Calf compresses (water approx. 30 - 35°C, renew every 15 minutes).

Important: Do not lower the temperature below 38°C, drink a lot and rest in bed.

Should these measures not reduce the temperature, then a doctor should definitely be consulted.

CLEANSING AND DISINFECTION

This thermometer is guaranteed waterproof and can be dipped into water or disinfection liquid. For disinfection please use one of the below recommended substances.

Descosept AF, Dr. Schumacher GmbH
Quick disinfection for wiping and spraying

Lyselol FF Concentrate, S+M
(dilution with water)

Ethyl alcohol 95.0%
(dilution with water); dip in for 24 minutes.

Advice: Any other cleansing substances, dilutions or cleansing method can cause disfunction!

CHANGING THE BATTERY

As soon as the symbol "■" appears in the LCD, the battery must be replaced.

- Remove the cover of the battery compartment (cap).
- Loosen the battery holder.
- Insert new battery (type LR 41 Alkaline, 1.55 V) (ensuring that it is the right way around, with the "+" facing upwards).
- Pay attention that by putting the lid back onto the battery compartment you won't damage the casing. Batteries are obtainable from all dealers specializing in electrical appliances.

Caution: Empty batteries and technical devices do not belong in the domestic refuse. Please dispose them at a recycle point in your area.

Technical data

Typ: Maximum thermometer
Measuring range: 32.00°C to 44.00°C
Temperature under 32.00°C
display "Lo" for low (too low)
Temperature over 44.00°C
display "Hi" for high (too high)

Accuracy: ± 0.1°C
between 35.50°C and 42.00°C
with surrounding temperature between 15°C and 25°C

Display: Liquid crystal display (LCD) with 4 digits, smallest display unit 0.01 °C

Acoustic Signal: when temperature increase less than 0.1 °C within 16 sec

Storage temperature: -10°C to +60°C

Battery: Type LR 41, 1.55 V

Battery check: Optic display ("■")

UEBE Medical GmbH
97877 Wertheim/Germany

CE 0123

• Directive for Medical Devices (93/42/EWG)
• DIN EN 12470-3: 2000 Compact electrical thermometers with maximum temperature indicator
• Protection grade against electric shock: Type BF

Note for hospitals, doctor's surgeries etc.:

In the event of professional use the measurements have to be checked every two years. This may be conducted by UEBE Medical GmbH or by an official agency responsible for calibration; a fee is charged for calibration.

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Manufacturer's guarantee

The equipment was manufactured and tested with great care. However, in the unlikely event of a defect being detected after delivery, we are prepared to give a guarantee in accordance with the following conditions:

1. During the guarantee period of 2 years from the date of purchase, we reserve the right to either repair any such defect at our expense or supply a perfect replacement unit (in any case the defective equipment must be returned to our factory).

2. Excluded from the guarantee are parts subject to normal wear and tear as well as damage caused by non-compliance with the operating instructions, improper handling (for example: unsuitable supply mains, breakages, leaking batteries) and/or dismantling of the unit by the purchaser. Furthermore, no claims for damages against us can be justified through the guarantee.

3. Claims under the guarantee can only be made during the guarantee period by presentation of the purchase receipt. In the case of claim under the guarantee, the equipment together with the purchase receipt and a written justification of the claim are to be sent to UEBE Medical GmbH, Service-Center, Zum Schlag 18, D-97877 Wertheim.

4. In the case of defectiveness of the goods, the contractual rights of the purchaser to claim against the supplier/seller in accordance with § 437 BGB (The German Civil Code) are not limited by the guarantee.

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